

CAMERA & SENSORS

HRI detects patient presence and supports interaction.

TOUCHSCREEN INTERFACE

HRI has simple, guided steps to enter symptoms and confirm details.

CONNECTED VITALS STATION

HRI contains a vital station that captures key vital signs through integrated/connected devices.

PROXIMITY SENSOR

HRI contains a built-in proximity sensor to detect when a patient approaches and wakes the system.



STATUS INDICATOR

HRI has a status light to indicate that the patient is finished and it's ready for the next patient.

HRI WORKFLOW IN MERCER

Patient arrives → approaches or joins the CARISURG pre-triage queue → robot detects and greets patient → language/accessibility needs selected → rapid safety screen → robot guides patient to the assessment area → symptom intake and vital signs collected → AI performs preliminary risk screening → information is sent to the nurse dashboard → nurse makes the final triage decision → patient receives queue/next-step guidance.

Urgent bypass: If red-flag symptoms are detected at any stage, the robot immediately alerts clinical staff and routine pre-triage is stopped.